

**BioMEDanas**  
FOR BIOMEDICAL ENGINEERING





BioMEDanas UG (haftungsbeschränkt)  
Fanny-Lewald-Ring 165, 21035 Hamburg  
Local court Hamburg  
Commercial register HRB 171489



Bank details  
Postbank NDL Deutsche Bank  
BIC: PBNKDEFF  
IBAN: DE38 1001 0010 0709 9601 34

CEO  
Eng. Anas Musleh  
tax number : 44/708/01234  
VAT ID No. DE347062075



# SK-P009

Double Arm Electric Surgical pendant



**Use:** ICU room, Suction ceiling installation.

**Material:** Main body use aluminium alloy material.  
Surface use powder-coated treatment.



## SK-P001

ICU Pedant Ceiling-Mounted Rail System(Cantilever)

Configurations:		SK-P001
Gas outlets standard	Oxygen	1 pc
	Vacuum	1 pc
	Medical Air	1 pc
Lighting lamps		2 pcs
Swing-arm lamp		1 pc
Electric sockets		6 pcs
Ground terminal		2 pcs
Net connect		RJ45 1 pc
Trays		3 pcs
Drawer		1 pc

Parameters:		SK-P001
Working Power		Customerization
Bridge length		2500-3200mm
Left and right move distance		300-600mm
Horizontal rotation angle		≤340°
Rotation angle		0-340°
Loading capacity		≤150kg
Gas outlets life		Insert&pull out more than 20000times

**Use:** ICU room, Suction ceiling installation.

**Material:** Main body use aluminium alloy material.  
Surface use powder-coated treatment.



**SK-P003**

ICU Pedant Ceiling-Mounted Rail System(Fixed Rotary)



**SK-P002**

Configurations:		SK-P002/SK-P003
Gas outlets standard	Oxygen	1 pc
	Vacuum	1 pc
	Medical Air	1 pc
Lighting lamps		2 pcs
Swing-arm lamp		1 pc
Electric sockets		6 pcs
Ground terminal		2 pcs
Net connect		RJ45 1 pc
Trays		3 pcs
Drawer		1 pc
Adjustable stainless steel iv pole		
Syringe pump support		1 pc
Storage basket		1 pc

Parameters:	SK-P002/SK-P003
Working Power	Customization
Bridge length	2500-3200mm
Left and right move distance	300-600mm
Horizontal rotation angle	≤340°
Rotation angle	0-340°
Loading capacity	≤150kg
Gas outlets life	Insert&pull out more than 20000times

**Use:** ICU room, Suction ceiling installation.

**Material:** Main body use aluminium alloy material,  
Surface use powder-coated treatment.



**SK-P004**

ICU Pedant Ceiling-Mounted Rail System



**SK-P005**

Configurations:		SK-P004
Gas terminal	Oxygen	1 pc
	Vacuum	1 pc
	Medical Air	1 pc
Instrument Platform		4 layers
Swing-arm lamp		1 pc
Lighting Lamps		2 pcs
Ground terminal		2 pcs
Power outlet		6 pcs
Network Interface		RJ45 1 pc
Drawer		1 pc

Parameters:		SK-P004
Working Power		Customization
Bridge length		1800-2200mm
Arm length		400mm
Horizontal rotation angle		≤340°
Loading capacity		≤150kg
Gas outlets life		Insert&pull out more than 20000times

**Use:** Operating room, Suction ceiling installation.  
**Material:** Use aluminium alloy material  
 Surface use powder-coated treatment



**SK-P008**  
 Double Arm Electric Surgical pendant



**SK-P009**

Configurations:		SK-P008
Gas connecting	Oxygen	1 pc
	Vacuum	1 pc
	Medical Air	1 pc
Rotating horizontal arm		1 pc
Electrical outlet		6 pcs
Gas braking system		1 set
Ground terminal		1 set
Adjustable stainless steel IV pole		1 pc
Adjusted ss injection pump shelf		1 pc

Parameters:		SK-P008
Working Power		Customization
One pc		
With three layers instrument platform		
Horizontal rotation angle		≤340°
Loading capacity		≤150kg
Gas outlets life		Insert&pull out more than 20000times

**Use:** General use, Suction ceiling installation.  
**Material:** Use aluminium alloy material  
 Surface use powder-coated treatment



**SK-P010**  
 Medical Multi-function Column



**SK-P011**  
 Medical Floor Stand Column

Configurations:		SK-P010
Gas outlet standard	Oxygen	2 pcs
	Vacuum	2 pcs
	Medical Air	2 pcs
Electric sockets		12 pcs
Adjustable stainless steel iv pole		1 pc
Layers		3 pcs
Syringe pump support		1 set
Drawer		1 pc
Storage basket		2 pcs

Parameters:	SK-P010
Working Power	Customerization
Horizontal rotation angle	≤340°
Loading capacity	200KG
Round corner anti-collision design	
Gas outlets life	Insert&pull out more than 20000times



**Use:** Operating room, ceiling type.  
**Material:** Use aluminium alloy material  
 Surface use powder-coated treatment



**SK-P014**  
 Single Arm Surgical Pendant(Adjustable Height)

Configurations:		SK-P014
Gas interface	Oxygen	2 pcs
	Vacuum	2 pcs
Plug		8 pcs
Ground terminal		2 pcs
Net work port RJ45		1 set
Adjustable ss IV pole		1 pc
Instrument platform		2 layers
Syringe pump support		1 set
Drawer		1 pc



**SK-P013**  
 Single Arm Electric Surgical Pendant



**SK-P018**  
 Single Arm Surgical Pendant (Adjustable Height)

Parameters:	SK-P014
Working Power	Customization
Working platform vertical lift range	600mm
Horizontal arm moving range	≤1000mm
Horizontal rotation angle	≤340°
Gas outlets life	Insert&pull out more than 20000times

**Use:** Operating room, ceiling type.

**Material:** Use aluminium alloy material  
Surface use powder-coated treatment



### SK-P012

Double Arm Surgical Pendant



### SK-P016

Double Arm Electric Surgical pendant  
(height Adjustable)

Configurations:		SK-P012
Gas interface	Oxygen	2 pcs
	Vacuum	2 pcs
Electricity socket		8 pcs
Ground terminal		2 pcs
Instrument platform layers		2 pcs
Net work port RJ45		1 set
Adjustable ss iv pole		1 pc
Syringe pump support		1 set
Drawer		1 pc

Parameters:		SK-P012
Working Power		Customization
Workbench vertical lift		600mm
Radius of cross arm		1400mm
Horizontal rotation angle		≤340°
Gas outlets life		Insert&pull out more than 20000times

**Use:** Abdominal Cavity , Ceiling type.

**Material:** Use aluminium alloy material  
Surface use powder-coated treatment



### SK-P020

Single Arm Electric Abdominal Cavity Pendant (Adjustable Height)



### SK-P021

Double Arm Abdominal Cavity Pendant  
(Fixed Height)

Configurations:		SK-P020
Gas interface	Oxygen	2 pcs
	Vacuum	2 pcs
	Air	1 pc
	CO2	1 pc
Electrical outlet		8 pcs
Ground terminal		2 pcs
Internet connector RJ45		1 pc
Adjustable stainless steel iv pole		1 pc
Syringe pump support		1 set
Drawer		1 pc
Syringe pump support		1 set
Drawer		1 pc

Parameters:		SK-P020
Working Power		Customization
Platform vertical lift range		600mm
Horizontal arm moving range		≤1000mm
Horizontal rotation angle		≤340°
With gas braking system		
Loading capacity		1400mm
With two layers instrument platform		
Gas outlet insert and draw life		20000 times

**Use:** Anesthesia, Ceiling type.

**Material:** Use aluminium alloy material

Surface use powder-coated treatment



### SK-P023

Single Arm Electric Anesthesia Pendant (Electric Height Adjustable)

#### Configurations:

		SK-P023
Gas interface	Oxygen	2 pcs
	Vacuum	2 pcs
Electric Sockets		8 pcs
Ground terminal		2 pcs
Net work port RJ45		1 pc
Instrument platform		2 layers
Suspension anesthesia machine		
Adjustable s <sub>i</sub> iv pole		1 pc
Syringe pump support		1 set

#### Parameters:

		SK-P023
Working Power		Customization
Workbench vertical lift		600mm
Lifting speed		22mm/s
Radius of the cross arm		≤1000mm
Horizontal rotation angle		≤340°
Gas outlet Insert&Draw Frequency		20000 times

**Use:** Endoscope, Ceiling type.

**Material:** Use aluminium alloy material  
Surface use powder-coated treatment



## SK-P024

Rotary Endoscope Pendant (Electric Height Adjustable)

Configurations:		SK-P024	Parameters:		SK-P024
Gas interface	Oxygen	2 pcs	Working Power	Customerization	
	Vacuum	2 pcs	Horizontal rotation angle	≤340°	
Electrical outlet		8 pcs	With gas braking system		
With equipotential ground terminal		2 pcs	Loading capacity	100-200KGS	
Internet connector RJ45		1 pc	With three layers instrument platform		
Adjustable stainless steel iv pole		1 pc			
Syringe pump support		1 set			
Drawer		1 pc			

**Use:** Endoscope, Ceiling type.

**Material:** Use aluminium alloy material  
Surface use powder-coated treatment



## SK-P025

Rotary Endoscope Pendant (Fixed Height)

Configurations:		SK-P025	Parameters:		SK-P025
Gas terminal	Oxygen	2 pcs	Working power	AC220V, 50HZ	
	Vacuum	2 pcs	Workbench vertical range	300mm	
Instrument platform	2 layers		Electric power	1 kw	
Electric socket	8 pcs		Lifting speed	22 mm/s	
Grounding terminal	2 pcs		Horizontal rotation angle	≤340°	
Adjustable ss iv pole	1 pc		Loading capacity	100kg	
Syringe pump support	1 set		Gas outle Insert & Draw Frequency	20000 times	
Drawer	1 pc				
Anti-collision round corner design					

Use: ICU/Infusion room, Suction ceiling installation.



**SK-P007**  
Infusion Rails System


Parameter : (Price/Meter)  
Working Power: Customization  
Loading capacity:  $\leq 150\text{kg}$   
Different color outlet prevent connect wrong  
Insert and pull out more than 20000times  
Electrical socket: according to drawbridge length  
Infusion hook: according to drawbridge length



**SK-P026**  
Rotary Surgical Pendant (Electric Height Adjustable)



**SK-P027**  
Rotary Surgical Pendant (Fixed Height)





**THANKS FOR  
WATCHING**

