

**BioMEDanas**  
FOR BIOMEDICAL ENGINEERING





BioMEDanas UG (haftungsbeschränkt)  
Fanny-Lewald-Ring 165, 21035 Hamburg  
Local court Hamburg  
Commercial register HRB 171489



Bank details  
Postbank NDL Deutsche Bank  
BIC: PBNKDEFF  
IBAN: DE38 1001 0010 0709 9601 34

CEO  
Eng. Anas Musleh  
tax number : 44/708/01234  
VAT ID No. DE347062075



## Advantage of Saikang patient monitor



Thickness: 8.9cm  
(thinnest in the field)

Unique fashion  
model design,  
high identification.



Industrial grade TFT screen,  
working life upon 50,000h  
(general less than 20,000h)



Lamp flashing &  
Tones alarm



Interface hidden design:  
dust cover



Import from OMRON company,  
life upon 1,000,000 times  
(domestic less than 500,000 times).  
Waterproof



Industrial shuttle imported from  
japanese ALPS company,  
working life upon 1,000,000h  
(domestic less than 100,000h)

## EM Series Integrated healthcare monitor



Models: EM032 12 for 12"  
EM032 10 for 10"  
EM032 8 for 8"  
Bed side/Healthcare department/Clinical

Basic configuration:  
EGC, RESP, HR, SPO2, PR, NIBP, TEMP

Optional configuration:  
Touch screen, 12 lead EGC, WIFI&GPRS, 2-IBP,  
ETC02, 5AH battery, Central monitoring system,  
Wall standing&Trolley

### Feature

- 12.1"/10.4"/8.4" high resolution TFT, color screen
- Fanless design
- ST and arrhythmia analysis
- SPO2 low perfusion index up to 0.1%
- Three alarm lights well distinguish technical and physiological alarm
- Pluggable Li-ion battery, max 5 hours continuous working time
- Connected with CMS by wire and wireless(3G and WIFI)
- Multi-language



## Patient Monitoring(ICU/CCU/OR)



Models: EM033 12S for 12"  
EM033 10S for 10"  
EM033 8S for 8"

ICU/CCU/OR



## Advantage of Saikang patient monitor



## Feature

- 12.1"/10.4"/8.4" high resolution TFT, color screen
- Semi-Modular and fanless design
- Multi-language
- Various of operation modes: touch screen, knob, mouse, keyboard
- 3/5 lead ECG
- ST and arrhythmia analysis
- STO2 low perfusion index up to 0.1%
- Three alarm lights well distinguish technical and physiological alarm
- Modular printer with 3 waveform output
- Pluggable Li-ion battery, max 5 hours continuous working time
- Connected with CMS by wire and wireless(3G and WIFI)

**Basic Configuration:**  
ECG, RESP, HR, SPO2, PR, NIBP, TEMP

**Optional Configuration:**  
Touch Screen, Recorder, 12 lead ECG, ETO2,  
2-IBP, WIFI&GPRS, 5AH battery, Wall standing  
&Trolley, Central monitoring system

## Central Monitoring System

Single/dual screen monitoring  
Monitoring 64 beds  
Wire and wireless(3G and WIFI)CMS  
Support HL7, HIS





## Pad Monitor

First aid, clinic examination, ward round, healthcare center, nursing homes, family healthcare




**Standard configuration:**  
SpO2/PR, ECG/HR, NIBP,  
Infrared forehead temperature.

**Optional configuration:**  
Blood glucose, CHOL(cholesterol),  
UA(uric acid)

## Feature

Android system, easy operation as iPad  
Small size & light weight  
10.1" TFT color touch screen  
Support Blood glucose, Forehead TEMP,  
Cholesterol, Uric acid.  
Rechargeable lithium battery, 10 hours  
continuous working time 120 hours data  
memory







**THANKS FOR  
WATCHING**

