

BioMEDanas
FOR BIOMEDICAL ENGINEERING





BioMEDanas UG (haftungsbeschränkt)
Fanny-Lewald-Ring 165, 21035 Hamburg
Local court Hamburg
Commercial register HRB 171489



Bank details
Postbank NDL Deutsche Bank
BIC: PBNKDEFF
IBAN: DE38 1001 0010 0709 9601 34

CEO
Eng. Anas Musleh
tax number : 44/708/01234
VAT ID No. DE347062075

The versatile ultrasound machine for large animal scanning

STANDARD
2
YEAR WARRANTY

Technology Driven
2-Year Warranty



Specification

Battery life	Lithium Ion 4~4.5 hours continuous use	USB 2.0	For connection to PC & Printer, download images to PC
Viewing options	I-Scan® Video goggles	Video	Composite Video (PAL&NTSC)
Frequency	2.0MHz to 10.0 MHz	Size	40 mm × 250 mm × 120mm
Digital channels	32	Weight	Only 0.6 Kg with battery and TFT LCD.
Cine-Loop	420~450	Temperature range	-10 °C to 40 °C / 14 °F to 104 °F *
Display	5.8 inches Color Transflective TFT		

Standard Configuration:

- 1*VetSan V3 machine
- 1* Abdominal Convex Probe C350 of 3.5 / 5MHz - (Ewes, goats, sows)
- 1* Li-ion high-power-capacity battery (5 hours/battery)
- 1*Battery charger (AC110V~250V,50/60Hz)
- 1*Neck strap ,hand strap
- 1*Aluminum alloy carry case
- 1*User guide
- 24 months warranty for machine and 12 months warranty for probes



Accessories – professional your VetSan V3



C350– Convex probe



C520– Microconvex probe



LH75– Linear probe



leather protector
(easy-clean by running water)



RC40V– Rectal convex probe



Rectal linear probe



I.R.F™ Introducer



I-Scan® V2
Video Goggles

The versatile ultrasound machine for large animal scanning

STANDARD
2
YEAR WARRANTY

Technology Driven
2-Year Warranty



Applications

VetSan V3

Thanks to 6 different probes , I.F.R™ introducer and I-Scan® video glasses, VetSan V3 can perform all veterinarian applications and even more easily. VetSan V3, featuring premium image quality, considerate design for comfortable, fast and reliable detection of pregnancy and gynecological diagnosis in difficult field conditions on a daily basis.

Maximum comfort while you work

- Only weight 0.6 kg;
- Excellent imaging quality and ergonomic design;
- Provides versatile solutions with the greatest ease of use for all veterinarian applications – bovine, equine, swine, ovine, Lama and small ruminants;
- 5.8 inch WVGA LCD monitor, wide viewing angle, LCD monitor can turn off when use I-Scan® Video Goggle
- 4.5~5.5 hours Li-ion battery for operating;
- Support USB 2.0 for ultrasound station to transfer images and video to your computer, 450~550-frame cine loop insider;
- Software & Report for reproductive system, and measurement for distance, area, circumference, volume, angle, heart rate;
- Uses a 3rd generation innovative navigation, direct access keys and user buttons for fast, single handed menu selection and setting changes;

Take VetSan V3 with you anytime and anywhere



- 10 hours continues working time (with 2 batteries)
- Robust, hard-wired probe designed for everyday use in the difficult field conditions



Portability by a shoulder belt and hand belt giving you total freedom of movements

The ultrasound scanner that makes you a specialist



I-Scan® video glasses eases your job and gets you focused (optional)



I.F.R™ introducer engineered to lessen fatigue and stress by allowing arm-free scanning, use in a variety of reproductive diagnostic procedures. (optional, compatible with the Rectal convex and Rectal Linear transducers)



Special design rectal convex probe based on customer feedback



Leather sheath (easy-clean by running water), designed to be resistant to extremely harsh farm conditions;




Good things come in small packages

Evaluate the reproductive tract
Foetal structures
Foetal sex

Superior Imaging Quality





**THANKS FOR
WATCHING**

